

Ketorolac on Lux 5µm Cellulose-3 in PO**Column:** Lux® 5 µm Cellulose-3, LC Column 250 x 4.6 mm, Ea**Dimensions:** 250 x 4.6 mm ID**Order No:** 00G-4493-E0**Elution Type:** Isocratic**Eluent A:** MeOH / IPA (90:10) Formic acid (0.1%)**Gradient****Profile:****Flow Rate:** 1 mL/min**Col. Temp.:** ambient**Detection:** UV-Vis Abs.-Variable Wave.(UV) @ 220 nm (ambient)**Analyst Note:** SecurityGuard™ Guard Cartridge System extends column lifetime.

- SecurityGuard Cartridges, Lux Cellulose-3 4 x 3.0 mm, 10/Pk Part No.: AJ0-8622

- Holder Part No.: KJ0-4282

20367

mAU

200

100

**Products used in this application:**

1

ANALYTES:**1** Ketorolac

Retention Time: 8.65 min

2 Ketorolac

Retention Time: 11.68 min

