

Extraction and Analysis of Metronidazole in Facial Cleansers**Column:** Kinetex® 2.6 µm XB-C18 100 Å, LC Column 50 x 2.1 mm, Ea**Dimensions:** 50 x 2.1 mm ID**Order No:** 00B-4496-AN**Elution Type:** Gradient**Eluent A:** 0.1% Formic Acid**Eluent B:** MeOH


Gradient Profile:	Step No.	Time (min)	Pct A	Pct B
	1	0	95	5
	2	2	60	40
	3	3	60	40
	4	3.1	95	5
	5	5	95	5

Flow Rate: 300 µL/min**Col. Temp.:** 25 °C**Detection:** Mass Spectrometer (MS) @ amu (650 °C)**Detector Info:** SCIEX< SecurityGuardSM ULTRA Guard Cartridge System extends column lifetime.

- SecurityGuard ULTRA Cartridges, UHPLC C18 for 2.1mm ID Columns, 3/Pk Part No.: AJ0-8782

- Holder Part No.: AJ0-9000

20631

**ANALYTES:****1** Metronidazole

Retention Time: 1.94 min

2 13-C2 15-N-Metronidazole

Retention Time: 1.93 min

**Products used in this application:**

Sample Preparation Details

for HPLC Application ID No.: 20631

Extraction and Analysis of Metronidazole in Facial Cleansers

PRODUCT DESCRIPTION:

StrataTM-X-C 33 μ m Polymeric Strong Cation, 30 mg / 3 mL, Tubes , 50/Pk

Order No.: 8B-S029-TBJ

SOLID PHASE EXTRACTION (SPE) PROCEDURE:

Note: The solvent volumes shown below are for a 30 mg bed mass.

The solvent volumes will need to be adjusted for a smaller or larger bed mass.

Condition:

Load:

- 1) Dissolve 0.250g of Neutrogena Foaming Face Cleaner With Microclear Technology in 10mLs of 0.1N HCL
 - 2) Vortex to achieve homogenous mixture
 - 3) Centrifuge for 5 minutes at 5000g
 - 4) Load
- 3mL aliquot of sample

Wash:

Dry:

5 min under full vacuum

Elute:

Final Prep and Analysis:

Evaporate to dryness at 50 C. Reconstitute in 200 μ L of 10:90 MeOH:0.1% Formic Acid

Inject: 10 μ L on HPLC Mass Spectrometer (MS) @ amu (650°C)

ANALYTES:	Spiked Conc. (ng/mL)	Log P	pKa	% Rec	%RSC (n=0)
1 Metronidazole	100				
2 13-C2 15-N-Metronidazole	1000				

Note: This method is designed as a convenient starting point for further investigation and can be tailored to meet your extraction goals.
Call your local Phenomenex Representative for assistance in method development and optimization techniques.

©2025 Phenomenex Inc. All rights reserved.

For more information contact your Phenomenex Representative at info@phenomenex.com



Phenomenex products are available worldwide.

www.phenomenex.com.cn

info@phenomenex.com