


Four Cations at 50ppb in the Presence of 1ppm Hydraxine

Column: N/A
Dimensions: N/A
Order No: N/A
Elution Type: N/A
Eluent A: 4 mM Nitric acid in 30% methanol
Gradient

Step No.	Time (min)	Pct A
1	0	100

Profile:
Flow Rate: 2 mL/min
Col. Temp.: ambient
Detection: Conductivity @ nm (ambient)
Analyst Note: SecurityGuard™ Guard Cartridge System extends column lifetime.

- SecurityGuard Cartridges, RP-1 4 x 3.0mm, 10/Pk Part No.: AJ0-5809
- Holder Part No.: KJ0-4282



ANALYTES:

- 1 Lithium
- 2 Sodium
- 3 Ammonium
- 4 Potassium
- 5 Hydrazine

