

Isopropanol Purity on ZB-624 30x0.32x1.8

Column: Zebron[™] ZB-624, GC Cap. Column 30 m x 0.32 mm x 1.80 µm, Ea
Phase: 6% Cyanopropylphenyl 94% Dimethylpolysiloxane
Dimensions: 30 meters x 0.32 mm x 1.8 µm
Order No: 7HM-G005-31
Oven Profile: 40°C for 4 min to 60°C at 5°C/min to 150 °C at 25°C/min for 2 min

Carrier Gas: Constant Flow Helium, 1.5 mL/min


Injection: Split 15:1 1 µL @ 150°C

Detection: Mass Selective (MSD) (0°C)

Analyst Note: Ethanol at 0.5% concentration was used as the internal standard. RSD values for water were calculated at 1.8% with n = 3. Detector was turned off for isopropanol solvent peak at 3.3 min to 4.0 min to prevent detector



Products used in this application:

**ANALYTES:**

- 1 Water
- 2 Ethanol

