


Carboxy-THC by MSD on a ZB-Drug1**Column:** Zebtron™ ZB-Drug-1, GC Cap. Column 10 m x 0.18 mm x 0.18 µm, Ea**Phase:****Dimensions:** 10 meters x 0.18 mm x 0.18 µm**Order No:** 7CD-G023-08**Oven Profile:** 220°C to 340°C @ 20 °C/min**Carrier Gas:** Constant Flow Helium, 55 cm/sec**Injection:** Split 10:1 0.1 µL @ 280°C**Detection:** Mass Selective (MSD) (320°C)**Analyst Note:** Analytes are 1.5ppm derivatized with BSTFA
Sim ions:
(Quant ion*) D3-Carboxy-thc (*374,476,477): Carboxy-thc (*371,473,474)**Products used in this application:**

17801

**ANALYTES:**

- 1 Carboxy-THC-d3
- 2 Carboxy-THC

