

NIDA-5 Extracted from Human Urine at Cutoff - Kinetex 2.6µm C18 50 x 2.1mm ID by LCMSMS

Column: Kinetex® 2.6 µm C18 100 Å, LC Column 50 x 2.1 mm, Ea

Dimensions: 50 x 2.1 mm ID

Order No: 00B-4462-AN

Elution Type: Gradient

Eluent A: 0.1% Formic acid in H₂O

Eluent B: 0.1% Formic acid in ACN

Gradient Profile:	Step No.	Time (min)	Pct A	Pct B
	1	0	95	5
	2	1.7	5	95
	3	1.9	5	95
	4	1.91	95	5
	5	2.6	95	5

Flow Rate: 1.2 mL/min

Col. Temp.: ambient

Detection: Mass Spectrometer (MS) @ amu (ambient)

Detector Info: ESI conditions

positive mode

NEB: 10.00


CUR: 11.00

IS: 4500.00

TEM: 450.00

CAD: 11.00

EP: 10.00



Products used in this application:



NIDA-5 Extracted from Human Urine at Cutoff - Kinetex 2.6µm C18 50 x 2.1mm ID by LCMSMS

ANALYTES:

- 1** Morphine
- 2** Morphine-d3
- 3** Amphetamine
- 4** Codeine
- 5** Codeine-d6
- 6** Methamphetamine
- 7** MDA
- 8** D5-MDA D6-oxycodone
- 9** MDMA
- 10** MDMA-d5
- 11** MDEA
- 12** MDEA-d5
- 13** Benzoylecgonine
- 14** Benzoylecgonine-d3
- 15** PCP
- 16** PCP-d5
- 17** THC-COOH
- 18** THC-COOH-d3

