


Metoprolol on Synergi Polar 4 μ m by SFC**Column:** SynergiTM 4 μ m Polar-RP 80 Å, LC Column 250 x 4.6 mm, Ea**Dimensions:** 250 x 4.6 mm ID**Order No:** 00G-4336-E0**Elution Type:** Gradient**Eluent A:** Methanol with 0.1% DEA**Eluent B:** Carbon Dioxide

Gradient Profile:	Step No.	Time (min)	Pct A	Pct B
	1	0	10	90
	2	2	10	90
	3	12	40	60
	4	14	40	60
	5	15	10	90

**Flow Rate:** 3 mL/min**Col. Temp.:** 35 °C**Detection:** UV-Vis Abs.-Fixed Wave. (UV) @ 254 nm (35 °C)**Detector Info:** Thar Discovery Series, Super Pure SFC instrument provided courtesy of Thar Technologies**Analyst Note:** SecurityGuardTM Guard Cartridge System extends column lifetime.

- SecurityGuard Cartridges, Silica 4 x 3.0mm, 10/Pk Part No.: AJ0-4348

- Holder Part No.: KJ0-4282

**ANALYTES:****1** Metoprolol