

## Aeris Ig-G comparison widepore XB-C18

**Column:** Aeris™ 3.6 µm WIDEPORE XB-C18 200 Å, LC Column 150 x 4.6 mm, Ea

**Dimensions:** 150 x 4.6 mm ID

**Order No:** 00F-4482-E0

**Elution Type:** Gradient

**Eluent A:** 0.1% TFA

**Eluent B:** 0.085%TFA in AcCN



| Gradient Profile: | Step No. | Time (min) | Pct A | Pct B |
|-------------------|----------|------------|-------|-------|
|                   | 1        | 0          | 90    | 10    |
|                   | 2        | 1          | 80    | 20    |
|                   | 3        | 15         | 40    | 60    |
|                   | 4        | 15.01      | 90    | 10    |

**Flow Rate:** 1

**Col. Temp.:** 40 °C

**Detection:** UV-Vis Abs.-Variable Wave.(UV) @ 214 nm (ambient)

22747



Products used in this application:



### ANALYTES:

1 Ig-G1 Standard

