

## Intact Reversed Phase Analysis of Infliximab-Temperature Effect (Detailed Overlay)

**Column:** bioZen™ 2.6 µm WidePore C4, LC Column 100 x 2.1 mm, Ea

**Dimensions:** 100 x 2.1 mm ID

**Order No:** 00D-4786-AN

**Elution Type:** Gradient

**Eluent A:** 0.1% TFA in Water


**Eluent B:** 0.1% TFA in Acetonitrile

| Gradient Profile: | Step No. | Time (min) | Pct A | Pct B |
|-------------------|----------|------------|-------|-------|
|                   | 1        | 0          | 80    | 20    |
|                   | 2        | 6          | 45    | 55    |
|                   | 3        | 6.1        | 10    | 90    |
|                   | 4        | 6.3        | 80    | 20    |
|                   | 5        | 6.5        | 10    | 90    |
|                   | 6        | 6.7        | 80    | 20    |
|                   | 7        | 6.9        | 10    | 90    |
|                   | 8        | 7.1        | 80    | 20    |
|                   | 9        | 8          | 80    | 20    |

**Flow Rate:** 0.5 mL/min

**Col. Temp.:** 70 °C

**Detection:** LC/UV (DAD, PDA) @ 280.000000000 nm (nanometers) (ambient)



## ANALYTES: